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EN

## **Part 1. Technical passport**


1. General description.....	24
2. Safety rules.....	25
3. Technical specifications.....	28
4. Package contents.....	28
5. Technical maintenance.....	30
6. Manufacturer's warranty.....	31

## **Part 2. Usage instruction**

1. General provisions.....	34
2. Conditions for treatment.....	37
3. Course of treatment.....	38
<i>Warranty certificate</i> .....	61
<i>Acceptance certificate</i> .....	64

# PART 1

## Technical passport

## 1. GENERAL DESCRIPTION

EN

External therapeutic foot electrode is intended for impact on reflexogenic zones of foot plantar surface.

The electrode is intended for multiple application in hospital, outpatient and home conditions.

The electrode is intended for joint application only with devices of dynamic electrostimulations.

**Compliance with standards:** *This medical device is CE marked according to the Directive 93/42 /EEC on medical equipment.*

## 2. SAFETY RULES

Please read the following safety instructions! They are designed to prevent user and property of damage. The manual contains important recommendations on proper use and care of the unit.

Patients with implanted electronic devices (such as heart pacemaker) and patients with idiosyncrasy to electric current are not allowed to use this unit.

The unit may not be used in the area of direct projection of the heart at the front.

Warning: Do not operate the unit while using another electrical medical device. This can cause burns and irreparable damage to the device.

Before you begin stimulation, make sure the patient is not connected to any high-frequency electronic device!

You must not switch the electrode to any other devices except the electrostimulators of DENAS and DiaDENS.

You must not switch two and more electrodes simultaneously to the one device.

**Attention!** If the patient is undergoing another course of physiotherapeutic treatment, use of this device is possible only after consultation with the attending physician.

The electrostimulation module must not be damaged.


The product contains fragile elements. Prevent it from shocks.


Avoid long-term exposure to direct sunlight at high (>25 °C) air temperature.


All work for repair of the product must be performed by qualified specialists at the company-manufacturer.

**Transportation conditions:** temperature from -50 to +50 °C, relative humidity from 30 to 93 %, atmospheric pressure from 70 to 106 kPa (from 525 to 795 mmHg).

**Storage conditions:** temperature from -50 to +40 °C, relative humidity from 30 to 93 %, atmospheric pressure from 70 to 106 kPa (from 525 to 795 mmHg).

**Operating conditions** temperature from +10 to +35 °C, relative humidity from 30 to 93 %, atmospheric pressure from 70 to 106 kPa (from 525 to 795 mmHg).

**Attention!** If the unit was stored at the ambient temperature below 10 °C keep it under normal climatic conditions for at least two hours before use.

**Recycling:** the unit packaging materials are not hazardous for the environment; they can be recycled.

The unit contains precious materials that may be used again after disposal with consideration of the environmental protection regulations.

They must be submitted to specially designated collection and recycling facilities (please consult in the relevant services in your locality).

**Attention!** *Any of controlled trials of DENS-therapy have not been carried out in pregnant women, so the device should be used **only in** case of absolute necessity, after consultation with DENS-therapy specialist.*

### 3. TECHNICAL SPECIFICATIONS

Specifications of the external therapeutic foot electrode are the following:


- weight of the device is no more than 500 g;
- external dimensions are no more than 370x370x50 mm;
- length of the electrode cable is not less than 50 mm.

### 4. PACKAGE CONTENTS

In table 1 you can see a list of delivery for the device.


*Table 1*

<b>Description</b>	<b>Quantity (pieces)</b>
External therapeutic electrode	1
Connecting cable	1
Consumer's container	1
Operation manual	1


**Fig. 1.** External therapeutic foot electrode DENAS-Reflexo.  
External view

Plug for con-  
nection to  
the electrode

Plug for con-  
nection to  
the de-  
vice

**Fig. 2.** Connecting cable intended for connection of the  
electrode to the devices of DENS-therapy:  
a) connection of the cable to the electrode DENAS-Reflexo;  
b) connection of the cable to the device

## 5. TECHNICAL MAINTENANCE

EN

Everyday technical maintenance includes the following procedures:

- external examination;
- disinfection (use standard disinfectants, soft lint-free and/or a small brush to clean the electrode);
- checking whether the electrode functions when connected to the devices of the series DENAS and DiaDENS.

## 6. MANUFACTURER'S WARRANTY

**6.1.** Company-manufacturer guarantees compliance of the device with technical conditions under observance of operating conditions, methods of transportation and storage.

**6.2.** Operating time of the device is 5 years. Under observance of operation regulations, actual operating time can significantly exceed the official.

**6.3.** Warranty period for the device is 24 months since the date of sale.

**6.4.** Completeness and appearance of the device is checked by the buyer in the seller's presence at the receipt of the product. Claims on completeness and appearance after the sale are not accepted.

**6.5.** In case of identification of defects during warranty period, the seller, manufacturer or manufacturer's representative specified by consumer rights legislation. The seller, manufacturer or manufacturer's representative is not responsible for the defaults should they occur after sale of the device in the result of:

- 1) failure on the part of the consumer to comply with the rules of transportation, storage, care and operation, provided by the present manual;
- 2) actions of third party;
- 3) force-majeure circumstances.

**6.6.** In case of unit defects or malfunction within the warranty period, the owner should send the following documents to the company-manufacturer's address: claim for the repair (replacement) with the name, address, telephone number; brief description of the malfunction, date and conditions of its appearance. All the defects or malfunctions are eliminated by the company-manufacturer or in service centers of the company-manufacturer.

***Address of the company-manufacturer:*****TRONITEK LLC**

620146, Russia, Ekaterinburg,  
Akademika Postovskogo Str., 15.  
Telephone: +7 (343) 267-23-30  
<http://www.denascorp.ru>  
e-mail: [corp@denascorp.ru](mailto:corp@denascorp.ru)


***Official representative in EU countries:*****DENAS-CZ s.r.o.**

VŘÍDELNÍ SŮL, spol. s r. o.  
Na Výhledě 886/3a, 36017 Karlovy Vary  
Czech Republic  
Tel./Fax: +420 353 549 285  
E-Mail: [denas@denasworld.com](mailto:denas@denasworld.com)  
[www.denasworld.com](http://www.denasworld.com)

# PART 2

## Usage instruction

## 1. GENERAL PROVISIONS

It is recommended to apply the external therapeutic electrode as an alternative to the inbuilt electrodes of the device for impact on reflexogenic zones of the soles. Application of the electrodes allows:

- to conduct treatment impact by yourself;
- to provide comfortable and convenient procedures at continuous and multiple treatment sessions;
- to make the conduct of the procedures easier.

Application areas: impact on feet is recommended at malfunction of lower limbs and also at diseases of internal parts of the body for stimulation of reflexogenic zones (fig. 3).


More detailed information on application of the device and also restrictions and contraindications are specified in the corresponding parts of the operation manuals of the devices for dynamic electrostimulation of nerves of the series DENAS and DiaDENS.

**Attention!** *It is prohibited to connect the electrode to any other devices, except devices of the series DENAS and DiaDENS.*

**Attention!** *Sudden pain in any location can be the first and sometimes the only sign of a serious disease. That is why if pain attacks appear for the first time and then repeat, and the intensity of the pain increases, it is necessary to consult a doctor urgently.*

**Fig. 3.** Projection of internal parts of the body to the foot plantar surface in the Su-Jok correspondence system

**Attention!** *In the case of breakdown of health during the treatment, you should consult with a doctor.*

**Attention!** *It is prohibited to use the electrodes in presence of solution of continuity of the skin (wounds, skin diseases) in the zone of impact.*

**Attention!** *Before using for the first time treat the electrodes of the device with a disinfection solution (e.g. 3% solution of hydrogen peroxide).*

**Attention!** *Lack of adequate cleaning of the elements can lead to allergic reactions or cross colonization - if the device is used by several people.*

**Attention!** *In case of appearance of any allergic reactions in the zone of contact of the device with the skin, you should stop using the device and consult a doctor immediately.*

**Attention!** *In case perceptibility of feet is very low, it is recommended to set the power in the range from 25% to maximum level.*

## 2. CONDITIONS FOR TREATMENT

EN

No special conditions are required for conducting procedures with the use of the external electrode. The room for the procedures should be dry and warm. During the session the patient sits in a position comfortable for him which ensures tight contact of the electrode with the skin.

### 3. COURSE OF TREATMENT

20 - 30 minutes before the session for hygienic purpose it is recommended to treat the work surface of the electrode with a soft brush or lint-free tissue, damped with a disinfection solution (e.g. 3% solution of hydrogen peroxide).

**Attention!** *If the skin is too dry it is allowed to damp slightly the skin surface with physiological solution 5 - 10 before the application of the external electrode to ensure tight contact and intensify treatment impact.*

**3.1.** Place the electrode on a flat dry surface (e.g. on the floor before the armchair).

**3.2.** Connect the electrode to the socket of the dynamic electrostimulation device (if necessary use a special adapter).

**3.3.** Put bare feet on the electrode so that the smooth dividing line appears to be approximately in the middle of the foot (fig. 4). Moving the feet on the surface of the electrode, find a position where the contact with the electrode is the tightest.


**Attention!** *The correct position of feet on the electrode is shown on figure 4.*

**Attention!** *Do not move the electrode during the stimulation!*

It is required to clean the contact part of the device before and after its use. Disinfection should be made with a standard disinfection solution (e.g. 3% solution of hydrogen peroxide).

**Fig. 4.** Correct position of feet on the electrode

**3.4.** Turn the device for dynamic electrostimulation on. Using buttons on the device set the necessary stimulation mode and power (intensity) of impact. The order for setting the mode and power is specified in the corresponding parts of the operation manual of the device for dynamic electrostimulation.

**3.5.** Conduct the session. The period of impact, necessity to treat additional zones, repetition of the procedures are defined individually and depend on the age, purposes and tasks of treatment. The recommended time of impact is specified in the operation manual of the device.

**3.6.** Upon termination of the session turn the device off. Disconnect the electrode and clean the work surface of the device (e.g. with 3% solution of hydrogen peroxide).

**Attention!** Keep the device in dry conditions.

A series of horizontal dotted lines for writing, spanning the width of the page.

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